


power carrousel™



PIPETMAN® M exclusive
charging stand



- | Rotating stand for up to five PIPETMAN M
- | Keep your pipettes organized
- | Fast charging saves on downtime



The POWER CARROUSEL is a comfortable charging solution for up to five PIPETMAN M (all models) that provides optimum protection from shocks and contamination and enables uninterrupted work!

Improves Your Efficiency

Save time with the POWER CARROUSEL and charge all of your PIPETMAN M at once in only one hour (to 80% of its capacity and 100% in three hours). You will not run out of battery!

Eliminate Clutter and Improve Organization

The slim design of the rotating POWER CARROUSEL conserves bench space for a more organized work area. Just rotate the carousel to reach the pipette that you need quickly and easily. The power cable is spooled and contained beneath the unit so you only need to expose as much cabling as you need.



MODEL	PART NUMBER
POWER CARROUSEL	FB1001

Gilson S.A.S.
19, avenue des Entrepreneurs BP 145 • F-95400
Villiers-le-Bel, France
T +33 (0) 1 34 29 50 00 • F +33 (0) 1 34 29 50 20

www.gilson.com/ContactUs

Gilson, Inc.
3000 Parmenter Street • PO Box 620027
Middleton, WI 53562-0027, USA
T 608-836-1551 or 800-445-7661 • F 608-831-4451

ISO 9001
CERTIFIED

LT900513/B © 03/2017 Gilson SAS.
All rights reserved. All trademarks
are the property of Gilson Inc.
and its subsidiaries.

GILSON®