

AmpliPur®

EXPERT TIPS

STERILE FILTERED PIPETTE TIPS FOR USE IN BIOLOGICAL APPLICATIONS



ISO Compliant Sterility

Manufactured from 94%
Recyclable Materials¹

Excellent Fit With added
Convenience



Sterility and Cleanliness for Reliable Results

Ensure the integrity of your experiments and the reliability of your data. The AmpliPur® Expert Tips, sterile version are certified with a Sterility Assurance Level (SAL) of 10⁻⁶ and are free from Human DNA/RNA and DNase/RNase. Ideal for microbiology, molecular biology, biomedicine, cell culture, and biobanking, they provide the high level of cleanliness and sterility required for your most critical applications.

Sustainable Design for a Greener Future

AmpliPur Expert Tips and their packaging are made from 94% recyclable materials¹, supporting your commitment to sustainability while reducing environmental impact. This eco-friendly design ensures that your work not only meets the highest scientific standards but also contributes to a more sustainable future.



Excellent Fit With added Convenience

AmpliPur Expert Tips, sterile version are designed to fit most common pipettes* found in the laboratory (ISO 8655 compliant). We are now elevating the Quality Standard, as well as extending the convenience and usability of the AmpliPur Expert Tips to address the needs of advanced biological workflows.

| Technical Information | |
|---------------------------------|---|
| Gilson Quality System | ISO 9001 |
| Product Compliance | ISO 8655, ISO 11137 |
| Certified Standards/Cleanliness | <ul style="list-style-type: none"> E-beam sterilized providing a Sterility Assurance Level (SAL) of 10⁻⁶ Free of Human DNA/RNA Free of DNase/RNase Hazardous Materials (2015/863/EU) |
| Materials (Recyclable) | Tips: Polypropylene (PP) |
| | Filter: Polyethylene (PE) |
| | Tray: Polypropylene (PP) |
| | Rack: Polypropylene (PP) |
| | Rack lid: Polypropylene (PP) |
| Materials (Disposable) | Plastic pouch: (Polyamide (PA) / Polyethylene (PE)) |
| Commercial Unit Packaging | 10 racks of 96 tips (960 tips) |
| Sterility Guarantee | Up to 5 years (depending on the expiry date) |

AmpliPur® Expert Tips, Sterile Version

| Volume Range | Nominal Volume of Compatible Pipettes** | Tips Qty/Rack | Tips Qty/ Commercial Unit (10 Racks of 96 Tips) | Tips qty/Overpack | Commercial Unit Part Number |
|--------------|---|---------------|---|-------------------|-----------------------------|
| 0.5-10 µL | 2 µL / 10 µL | 96 | 960 | 4800 | F174100 |
| 1-20 µL | 20 µL | 96 | 960 | 4800 | F174200 |
| 10-200 µL | 100 µL*** / 200 µL / 300 µL | 96 | 960 | 4800 | F174300 |
| 100-1000 µL | 1000 µL | 96 | 960 | 4800 | F174400 |
| 100-1200 µL | 1000 µL / 1200 µL | 96 | 960 | 4800 | F174500 |

*Pipette compatibility needs to be verified by the user.

**Within ISO 8655 specifications.

***Compatible with many 100 µL pipettes, except Gilson P100 single channel.

The recyclability of the materials is considered only when the materials are used to manufacture the product. It does not consider any end-use scenarios, or any contamination resulting from it. Recycling of a used product should be in accordance with local health and safety protocols, as well as guidelines that consider the impact of the contaminants on the environment. For further information, please consult your health and safety manager, environmental/sustainability manager, or local authority responsible for recycling lab waste.

gilson.com/contactus

LT800623/B | ©2025 Gilson S.A.S.
All rights reserved. All trademarks are the property of Gilson, Inc. and its subsidiaries.

