

MyPipetman®

ENTERPRISE

A UNIQUE COMBINATION OF FUNCTIONAL CHOICES AND PIPETTE PERSONALIZATION

Air-Displacement Pipettes

Design Your Very Own Gilson Pipette Series

MyPIPETMAN® Enterprise combines reliability, ergonomics, technical, and personalization options to make it the ideal solution to ensure efficiency and answer your pipetting needs within various workflows.

- **Select your preferred technical features** to optimize the use of your pipette and overall productivity in the lab.
- **Color choices and printing options** on your pipette allow for better organization and identification within your workflows.

1 Choose Your Technical Features

Volume Locking

- 3-position locking system
- Simple ratchet

Pipetting Forces (Single Channel Models)

- Reduced
- Standard

Tip Ejection Forces

- Reduced
- Standard

Tip Ejector (single channel, P2 to P1000 models)

- Metal
- Plastic



2 Choose Your Colors

Push-button (Lower part only)



Volume Locking and Ejection System



Handle



Window



Ring



Lower Parts



3 Add Extra Personalization to Your Pipettes

A

Personalize one side of your pipette with your company name or the text of your choice





Easy Customization

MyPIPETMAN Enterprise pipettes are made-to-order based on your design and technical choices.


Create and receive your custom pipettes in a few simple steps:

- Define your customization options with your Gilson Account Manager
- Confirm your choices
- We will prepare your pipettes with great care

MyPIPETMAN Enterprise pipettes are produced in a series of a minimum of **150 units** (mix of different volumes). All pipettes inside a series will have the same technical features, same design, and same information printed on the handle.



Quality Tips & Service



PIPETMAN DIAMOND® Tips are designed to fit your MyPIPETMAN perfectly and guarantee maximum accuracy and precision.

Benefit from high-quality pipette service with our flexible, tiered service level offerings, tailored to fit your lab's needs while increasing the lifespan of your pipettes. We can provide ISO 17025 compliant service in our dedicated Service Centers or on-site in your lab, perfect for GLP/GMP requirements.

Accessories



Other accessories like colored windows, tip holders, and tip ejectors are available in different colors and material for additional customization and ease of use.

Standard accessories, such as pipette stands and holders, are also available.

Please visit www.gilson.com for ordering information

MyPIPETMAN® Enterprise				
Model	PIPETMAN® DIAMOND Tips		Part Number	Volume (µL)
Single Channel Models				
P2	D10 DL10	DF10ST DFL10ST	FP10001E	0.2 - 2
P10	D10 DL10	DF10ST DFL10T	FP10002E	1 - 10
P20	D200	DF30ST	FP10003E	2 - 20
P100	D200	DF100ST	FP10004E	10 - 100
P200	D200 D300	DF200ST DF300ST	FP10005E	20 - 200
P1000	D1000 D1200	DF1000ST DF1200ST	FP10006E	100 - 1000
P5000	D5000		FP10007E	500 - 5000
P10mL	D10mL		FP10008E	1000 - 10000
Multichannel Models				
P8x10	D10 DL10	DF10ST DFL10ST	FP10009E	0.5 - 10
P12x10			FP10010E	
P8x20	DL10 D200	DFL10ST DF30ST	FP10011E	2 - 20
P12x20			FP10012E	
P8x200	D200 D300	DF200ST DF300ST	FP10013E	20 - 200
P12x200			FP10014E	
P8x300	D200 D300	DF200ST DF300ST	FP10015E	20 - 300
P12x300			FP10016E	

gilson.com/contactus

LT800619/A | ©2022 Gilson S.A.S.
All rights reserved. All trademarks are the property of Gilson, Inc. and its subsidiaries.

